Montage- und Gebrauchsanweisung	04	DEU
Installation and operating instructions	14	ENG
Istruzioni per il montaggio e l'uso	24	FRA
Mode d'emploi	34	ITA
Bruksanvisning	44	SWE

D^{care}®

AMALIA 2

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Care luminaire

Lampe de soins

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SYMBOLE; SYMBOLS; SYMBOLES; SIMBOLI; SYMBOLER

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
IMPORTANT!
THIS OPERATING MANUAL MUST BE THOROUGHLY READ BEFORE USING THE PRODUCT!
→ STORE FOR LATER USE!

CONTENTS

1.	VARIANTS AND SCOPE OF DELIVERY	15
1.1	AMALIA 2 P S4	15
1.2	AMALIA 2 P S7	15
1.3	AMALIA 2 B S8	15
2.	SAFETY INSTRUCTIONS	16
2.1	Intended use	16
2.2	User profile	16
2.3	Safety instructions	16
2.4	Warning levels	16
3.	MOUNTING	16
3.1	Workload data	16
4.	OPERATION	17
5.	CLEANING	18
6.	SAFETY CONTROLS	18
7.	DEMOUNTING	18
7.1	Disposal	18
8.	ACCESSORIES	19
9.	ADDITIONAL INFORMATION	19
10.	TROUBLESHOOTING	19
11.	TECHNICAL DATA	20
12.	ELECTROMAGNETIC COMPATIBILITY (EMC)	21

1. VARIANTS AND SCOPE OF DELIVERY


1.1 AMALIA 2 P S4

A: Lamp with ankle joint, gooseneck and pin

B: Plug-in power supply

1.2 AMALIA 2 P S7

A: Lamp with a double gooseneck and pin

B: Plug-in power supply

1.3 AMALIA 2 B S8

A: Lamp with double gooseneck and bed adapter

B: Plug-in power supply

2. SAFETY INSTRUCTIONS

2.1 Intended use

AMALIA 2 is an examination light. It is intended to illuminate the patient's body to support local diagnosis or treatment. Discontinuation of diagnosis or treatment because of a light failure may be possible at any time without risk to the patient.

12.2 User profile

Medical specialist

All individuals who have completed medical training and work in the professional field they trained for.

Cleaning specialist

Trained in national and workplace associated hygiene regulations.

Electrician

Trained in electronics and electrical technology and knows the relevant standards and regulations.

Qualified specialist

Qualified due to his technical training, knowledge and experience and knowledge of the rules, to perform assembly/dismantling.

2.3 Safety instructions

- ▶ Operation by a medical specialist
- ▶ This manual is part of the product and must be stored and made available to all future users.
- ▶ All work on the light (including repairs) may only be carried out by a qualified electrician. Installation may only be carried out by a qualified specialist.
- ▶ The manufacturer or a service technician authorized by the manufacturer or a similarly qualified person may only replace the light source of this luminaire.
- ▶ The light must not be altered or manipulated. Only approved original parts may be used. Any use other than the intended use with genuine parts may lead to other technical values and life-threatening hazards.
- ▶ Operation in potentially explosive areas is prohibited. The power supply to the light represents a potential ignition source.
- ▶ The light may only be operated in dry and dust-free rooms.
- ▶ The luminaire should not be left on without supervision.
- ▶ A damaged light should not be used. Defective cable locations also represent a potential hazard. Do not place cables near heat sources or on sharp edges.
- ▶ Never place extra weight on the light head and the arm system.
- ▶ This product emits potentially hazardous optical radiation. Do not look directly into the light emitted by the lamp. Eye injuries may occur.
- ▶ When in operation, the light must never be covered with a cloth or similar item.
- ▶ The ventilation openings (if available) must always be kept free during operation!

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2.4 Warning levels

DANGER

Warnings of dangers which, if the measures are ignored, can lead to death or serious injuries.

WARNING

Warning of hazards that can result in **injury** if instructions are not followed.

CAUTION

Warning of hazards that may cause **material damage** if instructions are not followed.

3. MOUNTING


3.1 Workload data

Bending moment M_B	25 Nm
Poids vertical F_G	90 N


3.2 AMALIA 2 P S4 and AMALIA 2 P S7

- ▶ **Fastening material is not included in the scope of delivery.**
- ▶ The lamp is equipped with a plug-in pin. The lamp must be connected to it in an accessory mentioned in chapter 8.
- ▶ If the "wall holder" interface is used, the installation must be carried out by qualified personnel.
- ▶ The wall must ensure a firm grip.
- ▶ Only use fixing material suitable for the corresponding substrate.

3.3 AMALIA 2 B S8

- Screw the threaded rod **B** into the light arm **A** until it stops and counter it with the nut **C**.

- Remove the cone **A** and nut **D** from the threaded rod **B**.
- Push the lamp with threaded rod **B** and cone **A** into the bed-side mounted cylinder **E**.
- Mount the cone **C** and the nut **D** and align the lamp.
- Tighten the nut **D**.

4. OPERATION

DANGER

Danger to life through an electric shock

Do not plug in a damaged power cable.

If there is any sign of damage to the power cable, replace it with a new one immediately


The connection voltage and frequency must match the data on the type plate.

CAUTION

Damage to property caused by an incorrect connection voltage.


Destruction of or damage to the lamp.

The connection voltage and frequency must match the data on the type plate.


- Plug in the cable
- Connect the plug-in power supply unit to the network
- All LEDs in the cone of light must light up before performing every function test.

4.2 Operation

- The 180° position (parking position) is not allowed.
- Sharp bends, approx. 90° on the connecting piece, are not allowed.

- An example of a correct position

WARNING

Eye injury warning

- This product may emit dangerous radiation. Never look directly into the light cone.
- The radiation emitted by this product complies with the exposure limit values for reduction of any risk of photobiological hazards, defined on the basis of IEC 62471.

- Switch the lamp on or off by briefly pressing the button.
- Press and hold the button to dim the lamp.
- Press and hold the button for 30 seconds to turn the power switch start function on or off.

5 CLEANING

DANGER

Danger to life through an electric shock

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CAUTION

Damage caused by improper cleaning

- ▶ Only such agents may be used for cleaning that do not impair the functionality of the light.
- ▶ No cleaning agents containing solvents, chlorine or abrasives may be used for cleaning, as such agents can lead to cracking in plastic parts, among other things.
- ▶ The agents used must be approved for use with plastics such as PC, PMMA, PA and ABS.
- ▶ Damage to the light due to concentrated disinfectants.
- ▶ For concentration and exposure time, please note the information on the supplementary sheet for the agent used.
- ▶ Scratches from use of unsuitable cloths.

RECOMMENDED DISINFECTANTS

- ▶ Lysoformin®
- ▶ Dismozon®
- ▶ Hexaquant®plus
- ▶ Sagrotan® quick disinfectant cleaner

CAUTION

Dirt reduces the luminosity

- ▶ Keep the screen clean by cleaning it regularly.
- ▶ Wipe clean only.

CAUTION

In order to minimise the risk of disease transmission, the applicable occupational health and safety regulations and the requirements of the national responsible committees for hygiene and disinfection must be observed in addition to these operating instructions.

6. SAFETY CONTROLS

DANGER

Danger to life through an electric shock

- ▶ When performing any work on the device, disconnect it from the power supply or pull the connector out of the socket and secure the device from being switched on again.

ANNUALLY:

- ▶ Check the connecting cable for damage and replace if necessary.
- ▶ Check for paint damage/cracks on plastic parts.
- ▶ Check for deformation or damage to the support system.
- ▶ Check for loosening of parts.

CAUTION

- ▶ Maintenance and repairs may only be carried out by qualified electricians.
- ▶ The corresponding user profile is in Section 1 Safety instructions.

7. DEMOUNTING

DANGER

Danger to life through an electric shock

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7.1 Disposal

Do not put the light in the household waste. Dispose of the light at a disposal point in accordance with local regulations or give it to a dealer with the appropriate service. Cut through the cable directly at the housing.

The products listed above are over 95% recyclable. For a high percentage of the used materials to be either physically reused or used for energy after the end of this product's life cycle, the luminaires have been designed with recycling in mind. They do not contain any dangerous substances or substances that require monitoring.

8. ACCESSORIES

Table clamp (Art. No.: D14.228.000 - black)

Universal fastener (Art. No.: D13.430.000 - pure white)

Rail clamp (Art. No.: D13.269.000)

Wall holder (Art. No. D13.231.000 - pure white)

Roller stand (Art. No.: D14.022.000)

Tubular holder (Art. No.: D12.313.000)

Bedside table adapter (Art. No.: D14.152.000)

9. ADDITIONAL INFORMATION

The light itself is maintenance-free.

Additional documentation for this product can be requested from the manufacturer on request.

The use of this light does not create any risks that could affect other devices.

To save energy, the light should only be switched on when it is actually being used.

Any serious incident that has occurred with the product **must** be **reported to the manufacturer or their representative and the responsible authorities of the country in which the user is located.**

The optical radiation emitted by this product complies with the exposure limits to reduce the risk of photobiological hazards according to IEC 60601-2-41.

10. TROUBLESHOOTING

Fault	Possible cause	Troubleshooting	User profile
Luminaire not working	Contact fault	Switch on again	All
Luminaire not working	No mains voltage	Check voltage, inspect all connections	Electrician
Luminaire not working	LED module defective	Contact the manufacturer	Only to be performed by the manufacturer's service personnel
Luminaire will not stay in position	Foot joint does not have enough friction	Adjust the ankle joint (1.25 Nm) until the desired friction is obtained.	Qualified specialist

11. TECHNICAL DATA

Electrical values:	
Rated input voltage	100 – 240 V
Frequency range	50 / 60 Hz
Power consumption at 100 – 240 VAC	2.3 - 2.5 W (5.7 - 8.9 VA)
Photometric values:	
Central illuminance E_c at a 1.0 m distance	1,200 lx *
Reference Distance, (D_{ref})	0.5 m
Light field diameter d_{10} at a distance of 0.5m	$\varnothing = 45 \text{ cm}^*$
Colour temperature	3'000 K *
Colour rendering index R_a	96 *
Colour reproduction index R_9	90 *
Ambient conditions for transport, storage and operation:	
Ambient temperature (storage and transport)	-20°C to +70°C
Ambient temperature (operation)	-10 °C hasta +35 °C
Rel. humidity (non-condensing) (storage and transport)	max. 90%
Rel. humidity (non-condensing) (operation)	max. 75%
Maximum utilisation height (during operation)	2000 m (above sea level)
Air pressure	80kPa
Maximum height for storage and transport	unlimited
Weight:	
AMALIA 2 P S4	1.5 kg
AMALIA 2 P S7	1.3 kg
AMALIA 2 B S8	1.7 kg
Operating mode:	
Type of operation	Continuous operation
Classification:	
AMALIA 2	Protection class II
Classification according to EU REGULATION 2017/745 (MDR), Item 51	Class I
U.S. FDA Device Class	Class I
Electrical safety testing and EMC according to:	AAMI ES60601-1: 2005/A2:2021 CAN/CSA-C22.2 No. 60601-1:14 IEC 60601-1:2005 + A1:2020 IEC 60601-1-2:2014 + A2:2020 IEC 60601-2-41:2021
Service life of the light source:	
Service life of the LED	60,000h (L70B50)

12. ELECTROMAGNETIC COMPATIBILITY (EMC)

Electrical medical devices are subject to special precautionary measures regarding electromagnetic compatibility. This device can be affected by other electrical devices.

WARNING

Use of this equipment in addition to or on other equipment should be avoided as this may result in improper operation. If such use is necessary, all devices involved should be observed to ensure they are functioning properly.

Use of accessories, transducers and cables not specified or provided by the manufacturer of this equipment may result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and may result in malfunction.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas and cables specified by the manufacturer) should be used no closer than 30 cm (12 inches) from any component of the AMALIA 2. Failure to do so may result in degradation of the performance of the device.

Electromagnetic environment

The device may only be operated in the environments which are defined in the section "Intended use" in the operating manual. The medical device is intended for use in the electromagnetic environment specified below

The models listed above are intended for use in the electromagnetic environment specified below. The customer or user of the models listed above should ensure that they are used in such an environment.		
Emission test	Compliance	Electromagnetic environment - instructions
RF emissions CISPR 11	Group 1	The models listed above use RF energy only for their internal operation. Therefore, their RF emissions are very low and are unlikely to cause interference with adjacent electronic equipment.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class C	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Compliant	

The models listed above are intended for use in the electromagnetic environment specified below.
The customer or user of the models listed above should ensure that they are used in such an environment.

Enclosure connector

Immunity test	Test condition	IEC 60601 Level of compliance	Electromagnetic environment - instructions
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact $\pm 2, 4, 8, 15$ kV air	± 8 kV contact $\pm 2, 4, 8, 15$ kV air	Floors should be made out of wood, concrete or ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Radiated RF EM fields and proximity fields from RF wireless communications equipment IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz: 80% AM 1 kHz	10 V/m 80 MHz to 2.7 GHz:	The mains power quality should be that of a professional healthcare facility and a home health environment.
	385 MHz (18 Hz pulse modulation)	27 V/m	
	450 MHz (FM+/-5 kHz deviation, 1 kHz sinusoidal or 18 Hz pulse modulation)	28 V/m	
	710 MHz (217 Hz PM)	9 V/m	
	745 MHz (217 Hz PM)	9 V/m	
	780 MHz (217 Hz PM)	9 V/m	
	810 MHz (18 Hz PM)	28 V/m	
	870 MHz (18 Hz PM)	28 V/m	
	930 MHz (18 Hz PM)	28 V/m	
	1720 MHz (217 Hz PM)	28 V/m	
	1845 MHz (217 Hz PM)	28 V/m	
	1970 MHz (217 Hz PM)	28 V/m	
	2450 MHz (217 Hz PM)	28 V/m	
	5240 MHz (217 Hz PM)	9 V/m	
	5500 MHz (217 Hz PM)	9 V/m	
	5785 MHz (217 Hz PM)	9 V/m	
RATED mains frequency magnetic fields IEC 61000-4-8	50 Hz or 60 Hz	30 A/m	Mains frequency magnetic fields should have the characteristics of a typical location in a typical commercial or clinical environment.
Immunity against nearby magnetic fields	30 kHz CW Test level: 8 A/m 134.2 kHz PM 2.1 kHz Test level: 65 A/m 13.56 MHz PM 50 kHz Test level: 7.5 A/m	Not applicable.	The AMALIA 2 lamp does not contain any magnetically sensitive components or circuits within the HOUSING.

The models listed above are intended for use in the electromagnetic environment specified below. The customer or user of the models listed above should ensure that they are used in such an environment.

AC CONNECTION input

Immunity test	Test condition	IEC 60601 Degree of Conformance	Electromagnetic environment - instructions
Fast transient electrical disturbances/bursts IEC 61000-4-4	± 2 kV 100 kHz repetition frequency	± 2 kV 100 kHz repetition frequency	The mains power quality should be that of a professional healthcare facility and a home health environment.
Surge voltage IEC 61000-4-5	± 1 kV voltage outer conductor - outer conductor ± 2 kV voltage outer conductor - earth	± 1 kV No applicable	The mains power quality should be that of a professional healthcare facility and a home health environment.
Conducted RF induced by RF fields IEC 61000-4-6	3 V 0.15 MHz – 80 MHz 6 V m) in ISM and amateur radio bands between 0.15 MHz and 80 MHz (n) 80% AM at 1 kHz	3 V 0.15 MHz – 80 MHz 6 V m) in ISM and amateur radio bands between 0.15 MHz and 80 MHz (n) 80% AM at 1 kHz	The mains power quality should be that of a professional healthcare facility and a home health environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% U_T ; 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0% U_T ; 0° 0% U_T ; 70% 0% U_T ; 0%	0.5 cycles 1 cycle 25/30 cycles (50/60 Hz) 250/300 cycles (50/60 Hz) (5 s)	The mains power quality should be that of a professional healthcare facility and a home health environment. If the user in the models listed above requires continuous force operation power outages, it is recommended that the models listed above be powered from an uninterruptible power supply or a battery power supply.

Note: n) The ISM (industrial, scientific and medical) - ranges between 0.15 MHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz. Amateur radio bands between 0.15 MHz and 80 MHz are 1.8 MHz to 2.0 MHz, 3.5 MHz to 4.0 MHz, 5.3 MHz to 5.4 MHz, 7 MHz to 7.3 MHz, 10.1 MHz to 10.15 MHz, 14 MHz to 14.2 MHz, 18.07 MHz to 18.17 MHz, 21.0 MHz to 21.4 MHz, 24.89 MHz to 24.99 MHz, 28.0 MHz to 29.7 MHz, and 50.0 MHz to 54.0 MHz.

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© Derungs Licht AG • D80.303.000 • 08.2025 • Index:1.0
Derungs Licht AG
Derungs Medical GmbH